



**Perfectly in tune**  
SATELEC® • Surgical • Tips

### BONE SURGERY

Osteotomy; Osteoplasty; Performing Bone Grafts.  
Allows cutting, excising and remodeling bone structures without any risk of soft tissue lesions.

<b>BS1</b> BS1 Slim- F87525 BL1 Long - F87527 BS1 - F87301	<b>BS2L</b> BS2L-2 - F87502 BS2L - F87302	<b>BS2R</b> BS2R-2 - F87503 BS2R - F87303	<b>BS4</b> BS4-2-F87504 BS4 - F87304	<b>BS5</b> BS5-2-F87505 BS5 - F87305	<b>BS6</b> BS6-2-F87506 BS6 - F87306
---	---	---	--	--	--

### SINUS LIFT

Lateral Sinus Lifts. A thin cut performed without effort, limiting any risk of sinus membrane perforation or antral artery damage. The elevation is performed by keeping a permanent contact between the membrane and the bone lining.

<b>SL1</b> SL1-2 - F87511 SL1 - F87311	<b>SL2</b> SL2-2 - F87512 SL2 - F87312	<b>SL3</b> SL3-2 - F87513 SL3 - F87313	<b>SL4</b> SL4-2 - F87514 SL4 - F87314	<b>SL5</b> SL5-2 - F87515 SL5 - F87315
--	--	--	--	--

### INTRALIFT


Sinus Lift by the Crestal Approach. The Intralift Kit makes it possible to perform hydronamically assisted sinus elevation with confidence and safety.

<b>TKW1</b> TKW1-2 - F87531 TKW1 - F87331	<b>TKW2</b> TKW2-2 - F87532 TKW2 - F87332	<b>TKW3</b> TKW3-2 - F87533 TKW3 - F87333	<b>TKW4</b> TKW4-2 - F87534 TKW4 - F87334	<b>TKW5</b> TKW5-2 - F87535 TKW5 - F87335
---	---	---	---	---

### EXTRACTION

Widening the Ligamentary Space. LC tips separate the tooth from its attachment system for a rapid atraumatic extraction, keeping the alveolar bone edge intact. Surgeries: extraction of wisdom teeth, fractured roots, impacted teeth, etc.

<b>Ninja</b> Ninja-2 - F87545 Ninja - F87345	<b>LC1</b> LC1-2 - F87507 LC1 - F87307	<b>LC1 90</b> LC1 90-2 - F87541 LC1 90 - F87341	<b>LC2</b> LC2-2 - F87542 LC2 - F87342	<b>LC2L</b> LC2L-2 - F87543 LC2L - F87343	<b>LC2R</b> LC2R-2 - F87544 LC2R - F87344
--	--	---	--	---	---



The most complete, efficient and secure range of Surgical Tips on the market.

**SATELEC®, World Leader**

- Inventor of the piezoelectric technology applied to dentistry
- Constant R&D investments in partnership with clinicians, universities and dental experts worldwide
- Products used daily in more than 100 countries

**SATELEC® tips, proven quality**

- Designed using high-tech CAD software
- Manufactured at the SATELEC® headquarters in Mérignac (Bordeaux, France) on high-precision equipment
- Strengthened by specific heat and surface treatment processes
- Individually tested to measure its efficiency, spray and acoustic performance

**Combine robustness and precision**

- Selection of high quality alloy for maximum resistance
- Control of each tip vibration in accordance with SATELEC® handpiece
- Longevity, robustness and reliability make the reputation of Piezotome®



### CROWN EXTENSION

Completely innocuous to the gingiva guaranteeing the integrity of bone and adjacent teeth.

<b>CE1</b> CE1-2 - F87551 CE1 - F87351	<b>CE2</b> CE2-2 - F87552 CE2 - F87352	<b>CE3</b> CE3-2 - F87553 CE3 - F87353	<b>BS6</b> BS6-2 - F87506 BS6 - F87306
--	--	--	--

### CREST SPLITTING

**NEW PRODUCT**  
The increasing thickness of the tips allows a smooth crest expansion limiting the risk of bone fracture.

<b>CS1</b> CS1-2 - F87561	<b>CS2</b> CS2-2 - F87562	<b>CS3</b> CS3-2 - F87563	<b>CS4</b> CS4-2 - F87564	<b>CS5</b> CS5-2 - F87565	<b>CS6</b> CS6-2 - F87566
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

**Surgical Tip Kit Promo**

Purchase any PiezotomeTip Kit\*  
(\*Bone Surgery (1&2) Kit, Sinus Lift (1&2) Kit, Intralift (1&2) Kit, Extraction (1&2) Kit or Crown Extension (1&2) Kits)

Get a BS1 Slim (Generation 1 or 2) Tip!  
BS1 (F87301) or BS1-2 Slim (F87525)

**\$253.00 Value**

**FREE!**

(See back cover for redemption rules and restrictions)