



**CROWN
EXTENSION**

ACTEON North America, SATELEC & SOPRO Companies Introduces the Crown Extension Kit

The objective of crown lengthening is to increase the clinical crown height of a tooth intended for a conservative or prosthetic restoration.

Compatible with Piezotome, Piezotome 2, and Implant Center2 piezoelectric bone surgery devices, the Crown Extension kit by SATELEC offers a complete solution to address the different stages of a surgical crown lengthening procedure. Specially developed for this procedure, the BS6 tip and three diamond-coated CE tips are designed for osteotomy (bone removal) and osteoplasty (bone reshaping).

The BS6 tip is used to perform larger shaping of the bone not supporting teeth. The CE1, CE2, and CE3 tips are all diamond coated spherical and cylindrical tips designed for very fine osteoplasty in interproximal, palatal, and buccal areas. Active on hard tissues and inactive on soft tissues, these tips offer a high level of safety when used close to delicate anatomical structures such as ligaments, nerves and arteries. This cutting selectivity renders the tips completely harmless to the gingiva, which is particularly beneficial following a gingivectomy with an internal bevel.

148 words

ACTEON North America
124 Gaither Drive
Suite 140
Mount Laurel, NJ 08054
(800) 289 6367 / (856) 222 4726
info@us.acteongroup.com
www.us.acteongroup.com

