

November 9, 2010

PRESS RELEASE:

ACTEON, North America releases new Perfect Margin Shoulder™ tips for the Newtron Line

ACTEON, North America (Mt laurel, NJ) has introduced a new set of ultrasonic tips for their Newtron® line of piezoelectric ultrasonic devices. These new tips are dubbed Perfect Margin SHOULDER tips. They have the same indications as the original Perfect Margin tips for inlay and onlay preparation but the ends of the tips have a flat edge instead of a chamfered end. These tips are designed for precise finishing of the gingival margin prior to impression taking as well as offer dentists the considerable advantages of selective abrasion or cutting using ultrasonics: the preservation of soft tissue and highly precise finishing of the prepared surface. The ultrasonic frequency of Newtron piezoelectric generators makes the instruments active only on hard tissues. There is no possible risk of bleeding at the sulcular level.

The pricing remains the same as the original Perfect Margin tips and the product reference numbers are as follows

- PMS1 - F02254
- PMS2 - F02255
- PMS3 - F02256
- Kit: F00736 (includes PMS1, PMS2, PMS3, and PM4 as well as an autoclavable stand and universal wrench)

Please contact your local dealer or ACTEON for more information