

So easy to use... It's *automatic!*

- Auto-access: automatic aperture of the slot device via motion detector
- Auto-detect: **PSPIX** automatically recognizes the size of the inserted plate
- Auto-scan: automatic readout after insertion of the Imaging Plate (avg 5 sec.)
- Auto-erase: the Imaging Plate is automatically erased after scanning
- Auto-eject: the Imaging Plate is automatically ejected in its collector for immediate reuse
- Auto-optimization: **PSPIX** adjusts the brightness and contrast of the image automatically - Image quality is always optimized
- Auto-sleep: **PSPIX** automatically switches to the energy-saving mode when it is not used for a certain period of time

ONE TWO

You shoot

You scan



Auto detect..



Auto scan..



Auto eject!



All the advantages of conventional film; flexibility and thinness... yet digital and reusable!

Size 0	Size 1	Size 2	Size 3
22 x 31 mm	24 x 40 mm	31 x 41 mm	27 x 54 mm


The thinness and flexibility of an imaging plate provides easy positioning in the mouth. The plates are wireless, hence easy to use and their surface is 100% active. They are available in four different sizes.



Onboard display



100 %
ACTIVE AREA



* Basic version is free with purchase of PSPIX.
** Premium version is available via upgrade.